

Confirmatory/Supplemental Testing

Serological Test Name	Confirmatory/ Supplemental Test Name	Testing Lab	Sample Requirements
HCV 2.0	Alternate Serology Test	Creative Testing Solutions	1.0 ml serum or plasma stored at 2-8 C for up to 7 days or frozen
HBsAg	HbsAg Neutralization	Creative Testing Solutions	1.0 ml serum or plasma stored at 2-8 C for up to 7 days or frozen
Anti-Hbcore	NA	NA	NA
Anti-HIV 1/2	IFA (Immunofluorescence Assay for HIV-1 Antibodies) <i>If IFA is negative or indeterminate reflex to HIV-2 EIA</i>	Creative Testing Solutions	0.5 ml serum or plasma stored at 2-8 C for up to 7 days or frozen
HTLV I/II	Anti-HTLV-I/II MP Western Blot	Creative Testing Solutions	.5 ml serum or plasma stored at 2-8 C for up to 7 days or frozen
PK TP (Syphilis)	RPR	TBC	1.0 ml serum stored at 2-8 C for up to 5 days or frozen
	Captia Syphilis-G (for PK TP positive, RPR negative donors)	Creative Testing Solutions	0.5 ml serum stored at 2-8 C for up to 5 days or frozen
Chagas (CGS)	Abbott ESA Chagas Assay	Creative Testing Solutions	1.0 ml serum or plasma stored at 2-8 C for up to 7 days or frozen
HIV, HCV, HBV alternate NAT	Procleix HIV-1/HCV/HBV NAT	Creative Testing Solutions	1.0 ml plasma. Store at 2- 8 C up to 7 days or frozen